

User manual Socks-jet

I. Preparation



Step 1

Place the Socks-jet on a flat, dry surface. The hand grip points towards you and the open tube side to the front. Wear rubber gloves to protect your stockings.



Step 2

Pull the stocking over the tube body and make sure that the toe of the stocking is facing forward.



Step 3

Pull the stocking further down over the Socks-jet until the heel knit appears at the top of the tube body.

II. Donning



Step 1

Place the Socks-jet in front of you on the floor and stabilize your position or sit down. Slide your foot into the open stocking.



Step 2

Slide your foot down until your heel is correctly positioned in the stocking. Now pull the stocking upwards until it is clear of the tube. Please check the comfort of your toes and heel again.



Step 3

Massage your stockings to their final position and smooth them with your hands towards your knee or thigh. Make sure you have a wrinkle-free fit.

III. Doffing



Step 1

Place your foot in the Socks-jet and sit down. Slide the stocking over the tube body from above and slip it downwards as far as possible, past the clamping jaws.



Step 2

Use the hand grips and press them inwards, so that the clamping jaws push against the tube.



Step 3

Tilt the Socks-jet and push it straight forward. Now slowly pull your foot out of the stocking.

IV. Passive use



With Assistance

Use the integrated hand grip for a simple application by nursing and support staff.